



MDXUSA[®]
Medical Equipment
Solutions



Accuit Sign 6

Multi-parameter Monitor

Features

- 8.4" color TFT display
- Suitable for adult, pediatric and neonatal patient
- Basic parameters (ECG, SpO₂, NIBP) in a durable case for bedside monitoring and transport monitoring
- Audible and visual alarms with adjustable alarm ranges
- Networkable with central monitoring system
- Powerful data management and storage capacity

NIBP

- Method: Oscillometry
- Operation modes: Manual/Automatic/STAT
- Measurement unit: mmHg/kPa selectable
- Measurement types: Systolic, Diastolic, Mean
- Over-pressure protection: Yes

ECG

- Lead Mode: 5-lead
- Input: RA; LA; RL; LL; V
- Sweep speed: 12.5mm/s, 25mm/s, 50mm/s
- Accuracy: ± 1 bpm or $\pm 1\%$, which is greater
- Protection: Withstand 4000VAC/50Hz voltage in isolation against electrosurgical and defibrillation
- S-T detection: Yes, Arrhythmia analysis: Yes
- Alarm: Yes, audible and visual alarm, alarm events recallable

SpO₂

- Measuring range: 0-100%
- Accuracy: ± 2 digit (70~100%);
0~69% unspecified

Standard Configuration

- ECG, NIBP, SpO₂

Options

- RESP, 2-TEMP, 2-IBP, EtCO₂, Network, Thermal Printer




Accuit Sign 6

Accessories



IBP Cable



IBP Module



IBP Transducer



EtCO₂ Sampling tube



Printer



EtCO₂ Module (Sidestream)

MDX Instruments USA, Inc

Bridgeport, Connecticut
USA

Email: info@mdxinstruments.com
Web: www.mdxinstruments.com

MDX reserves the right of product appearance, specifications and function modification and updating, without prior notification. (10-07-01) V1.0